

USER MANUAL

Ver. 2.1



SGM-STIM-2 CONSTANT CURRENT, MULTIMODE STIMULATOR FOR DIRECT CORTICAL ELECTRICAL STIMULATION

SGM-STIM-2

WARRANTY

SGM d.o.o. warrants that our products will be free from defects in material and workmanship under normal use and service for a period of twenty-four (24) months from date of shipment from our warehouse.

This warranty is void if the Product is subject to accident, misuse, neglect, improper handling, improper installation, improper repair or is modified or altered. The sole and exclusive obligation of SGM under this warranty is the repair or replacement of such defective or missing parts which are causing the malfunction. If SGM does not replace or repair such parts, the end-user's sole remedy against SGM shall be to obtain a refund of the price paid for such Products(s) as are proven to be defective. To obtain service under this warranty, the end user must bring the malfunction to the attention of SGM within the twenty-four (24) month period following the product shipment from SGM and no later than thirty (30) days after the occurrence of such malfunction, whichever occurs first.

SGM shall not be liable under this warranty unless SGM's examination of the products shall disclose, to its satisfaction, that the defects have not been caused by misuse, neglect, improper handling, improper installation, improper repair, modification or alteration.

THERE ARE NO OTHER WARRANTIES WHICH EXTEND BEYOND THE DESCRIPTION ON THE FACE HEREOF. ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING THE IMPLIED WARRANTIES OF MERCHANTABILITY, AND FITNESS FOR A PARTICULAR PURPOSE, ARE HEREBY EXCLUDED.

TABLE OF CONTENTS

GENERAL PRECAUTIONS	4
WARNINGS AND CONTRAINDICATIONS	5
1. TECHNICAL DESCRIPTION	6
1.1 INTRODUCTION	6
1.2 STIMULATOR CONTROLS	7
1.2.1 STIMULATION PARAMETERS SETUP	7
Multipulse mode	7
Penfield mode	7
Output	7
Current out	7
Frequency	7
Pulse width	7
Timer	7
1.2.2 BACK PANEL	9
1.3 TECHNICAL SPECIFICATIONS	9
1.3.1 STIMULUS OUTPUT	9
1.3.2 STIMULUS PULSE RATES	10
1.3.3 TRIGGER FACILITIES	10
1.3.4 OTHER	10
1.4 MAINS CONNECTION	10
1.5 ENVIRONMENTAL CONDITIONS	11
1.6 OUTPUT IMPEDANCES AND CURRENT	11
1.7 ELECTRODES INTENDED FOR USE 2	11
2. INSTRUCTION FOR USE	12
2.1 PRECAUTIONS WHEN USING THE SGM-STIM-2 STIMULATOR	12
2.2 POWERING THE STIMULATOR	3
2.3 DEFINITION OF STIMULUS PROPERTIES.....	14
2.4 OPERATING MODE SETUP.....	14
2.5 CONNECTING THE PATIENT.....	14
2.6 STIMULATION	14
2.7 END OF STIMULATION.....	14
2.8 STIMULATING ELECTRODES.....	14
2.9 OUTPUT LIMITS	15
2.10 DEVICE INSPECTION AND MAINTENANCE	15
2.11 CLEANING	15
3. ACCESSORIES.....	16
4. EMC TABLE	18
SGM-STIM-2 - EMC TABLE	18

GENERAL PRECAUTIONS

Operation with SGM-STIM-2 stimulator is allowed only to a highly qualified operator. Please make sure that all users of the SGM-STIM-2 system read these instructions and understand the conditions for safe use of the system. Failure to understand the information may lead to dangerous situations.

Carefully read this User Manual before using the SGM-STIM-2 stimulator.

HV switch must always be easily accessed and should not be covered with anything.

If there is a DC component on the device output – immediately turn off the device, and send it to repair.

Protective earth terminal must always be properly connected.

Only equipment that complies with IEC 60601 should be connected to the trigger input or output of this stimulator. To prevent unpredictable effects, only devices complying with safety standard IEC 60601 are allowed in the patient environment.

Electrical equipment for medical uses requires special EMC precaution and needs to be installed and serviced according to the EMC documentation of the device.

SGM-STIM-2 is inappropriate for use in oxygen rich environment.

Portable and mobile RF equipment can affect electrical equipment for medical use.

SGM-STIM-2 is intended for use only in moderate climates of the surgery room.

No modification of this equipment is allowed.

SGM-STIM-2 should be disposed of in accordance with the European directive 202/96/EC.

Ensure that during operation SGM-STIM-2 is positioned in such way to allow operating of the disconnection device.

The patient connections are driven from an isolated current source. The leakage values from the patient connections to the ground reference are within relevant limits for class BF equipment according to IEC 60601 standard.

Device is classified as TYPE IIb CLASS BF equipment.

SGM-STIM-2 does not display any system messages during operation.

SGM-STIM-2 output extender is defined as applied part (see section 3 - Accessories)

For all servicing requirements please contact:

SGM d.o.o.

Grge Novaka 22a

21000 Split, Croatia

Tel: ++385 (0) 21 779 916

Email: sales@sgm.hr

Web: www.sgm.hr

WARNINGS AND CONTRAINDICATIONS



- Do not use SGM-STIM-2 on patients with cardiac and brain pacemakers
- Avoid using SGM-STIM-2 on areas of the body in close proximity to any implantable electric device
- SGM-STIM-2 should not be used for stimulation in the vicinity of uterus during pregnancy
- SGM-STIM-2 may cause involuntary spasms of muscles of the mouth, tongue and lower jaw
- Do not remove device covers - high voltage is present within the unit.
- Device is not to be used for transthoracic surgeries.
- Do not use this device in explosive atmosphere or in the presence of flammable anesthetics.
- Avoid accidental contact between connected but unapplied APPLIED PARTS and patient or other conductive parts including those connected to protective earth.
- SGM-STIM-2 is not defibrillation-proof.
- Connection of a patient to HF surgical equipment and to evoked response equipment simultaneously may result in burns at the site of the electrical stimulator and possible damage to the electrical stimulator
- Operation in close proximity (for example 1 m) to shortwave or microwave therapy equipment or GSM devices may produce instability in electrical stimulator output.
- The stimulator is not suitable for use in domestic establishments and in establishments directly connected to a low voltage power supply network which supplies buildings used for domestic purposes.

This stimulator is intended for use by healthcare professionals only.

This stimulator may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as re-orienting or reallocating the stimulator or shielding the location.

Standard sequence of operation is described in sections 2.2 to 2.7 of this document.

WARNINGS AND CONCLUSIONS FROM

**CLINICAL EVALUATION REPORT According to MEDDEV 2.7.1.
For the Medical Device: SGM – STIM – 2 (CORTEXSTIM)**

8.0 CONCLUSION

- A. Direct cortical subcortical electric stimulation using muscle responses (muscle MEP) provides an effective means of monitoring and/or mapping motor pathways during surgery. Although complete optimal warning criteria have not yet been elucidated, simple criteria such as disappearance of the response are very useful guides.
- B. Care must be taken to prevent injury and to assure that anesthesia does not affect the muscle MEP responses.
- C. Significant changes in the muscle MEP during surgery bear a strong correlation with injury to the motor pathways.
- D. Patient's risk of injury using this stimulator is very low, while the benefit in preventing and detecting neurologic injury is very high.

1. TECHNICAL DESCRIPTION

1.1 INTRODUCTION

SGM-STIM-2 stimulator is designed to be used for direct electrical stimulation of cerebral cortex focal regions for motor, sensory and language functions. Constant current types of stimulus are used to deliver controllable amount of electrical charge, regardless of the impedance of tissue between electrodes.

The stimulator can be used in classic Penfield mode, or it can use the novel Multipulse method for intraoperative brain stimulation. Operation mode has to be selected before definition of stimulation parameters. For each method, parameters are defined independently and are preserved during a session.

Dual functionality enables:

- comparison of stimulation efficiency for these two methods of stimulation;
- use of a more appropriate method for specific intraoperative application (e.g. motor area, speech area, cognitive functions).

Penfield or multipulse stimulation method can be alternatively used.

For both types of stimulation, basic parameters are adjustable before, and can be changed during the stimulation process.

Stimulation probes that can be used with this stimulator are eight-pole strip electrodes or handheld probe.

Different trigger modes of operation can be used:

- auto triggered
- local manual trigger
- external trigger

Stimulator can be used as a stand-alone device with full local functionality, or can be connected and remotely controlled from IOM system or remote computer.

1.2 STIMULATOR CONTROLS

The stimulator user interface is organized on the front and rear device panel. All controls that are used before the patient is connected to the stimulator are on the rear panel. These are: mains connection elements, input connector for triggering the device from an external source, and output connector for triggering other equipment. Controls on the front panel are used for adjustments and stimulation control during patient stimulation.

1.2.1 STIMULATION PARAMETERS SETUP

LCD display and simple navigation keyboard are used for stimulation parameters setup.


Up/Down arrow key enables vertical navigation within the display to select stimulation parameter. During this navigation, a small yellow circle indicates the selected parameter.

Left/Right arrow key is used for value adjustment of the selected parameter. After switching POWER-ON and a short period of booting, starting menu enables selection of one of stimulation modes:

- multipulse stimulation using strip electrode
- multipulse stimulation with hand-held electrode
- Penfield stimulation using strip electrode
- Penfield stimulation with hand-held electrode

For any of these modes, default setup of stimulator parameters is presented on the LCD display.

MULTIPULSE MODE	
Output	1
Current out	0 mA
Pulses	1
Trigger	Auto 1.0 Hz
Interstimulus	4.0 ms
Pulse width	500 μ s
PENFIELD MODE	
Output	1
Current out	0 mA
Frequency	50 Hz
Pulse width	1000 μ s
Timer	0 s



Output selector: Second row on the presentation panel shows the selection of hand-held or selected active electrode in a strip montage. Using Left/Right key, active electrode can be selected.

Current out: Menu entry is used for output current intensity adjustment - it can be increased or decreased from 0 to 40 mA using right or left key. For output current from 0 to 5 mA one step increment is 0.2 mA. For higher current output, 1 mA increment step is used.

Pulses: Menu option enables selection from 1 to 9 output pulses within one train.

Trigger: “Auto” option is preset as default, with 1.0 Hz repetition rate. This mode is used for continuous stimulation with single pulse or group of pulses at repetition rate between 0.1 Hz and 8.0 Hz with 0.1 Hz resolution. Stimulation can be started and stopped with the yellow PULSE - START- STOP button.

Using left navigation arrow key, trigger mode can be adjusted to “external” or “manual” mode. “External” triggering mode will use trigger pulses from other electronic equipment to generate output from the stimulator. “Manual” triggering mode - allows operator to send single output stimulus (pulse or group) as determined by depressing PULSE-START-STOP switch on the front panel.

With similar manipulation, interstimulus and pulse width menu entries can be used for adjustment of interstimulus interval between 1 and 10 ms, and pulse width from 100 to 1000 μ s. To prevent “overlapping” of output pulses within one group, if pulse width is set to 1000 μ s, interstimulus interval cannot be decreased below 1.1 ms. Also, if interstimulus interval is adjusted to 1.0 ms, maximum pulse width is limited to 900 μ s.

OUTPUT CONTROL

This section is used for controlling the delivery of output pulses to the patient. It consists of HIGH VOLTAGE ON/OFF switch, PULSE START/STOP push button and multipole connector for strip or hand-held electrodes connection.

Black toggle switch in the HIGH VOLTAGE ON position enables the unit to deliver output pulses. This is not the primary power supply switch; it blocks functioning of the high voltage module. HIGH VOLTAGE OFF position disables output pulses while the trigger output will still deliver output pulses.

OUTPUT INDICATION:

Green indicator in the HIGH VOLTAGE ON/OFF switch shows that the internal high-voltage supply is activated, and output stimulus will be delivered. If output current cannot be delivered (due to wrong output selection, disconnection of electrode wiring or disconnection of electrode from the patient), this indicator will start to flash at high frequency.

PULSE START/STOP push button is the manual trigger switch, or Foot Pedal. Its function depends on the selected trigger mode:


In manual trigger mode, upon depressing this switch, a single output pulse or group of pulses would be produced, characterized by the pulse settings.

In auto trigger mode, depressing this switch will activate or stop automatic generation of output pulses. After each activation of this control, red light indicator (on PULSE START/STOP switch) shows that stimulus is sent to the patient.

In external trigger mode, PULSE START/STOP push button is disabled and delivering of pulses is controlled by external triggering device.

1.2.2 BACK PANEL

On the back panel, mains connection and triggering sockets are located. Mains inlet, power switch and mains fusing are integrated in one block. Potential equalizing DIN 42801 stud is used as a potential reference for unit and bonding point.



TRIGGER OUT socket delivers positive TTL pulse in synchronism with the start of each stimulation pulse in a single pulse mode of operation. When pulse groups are being generated, one output trigger pulse is produced in synchronization to the first pulse in the train.



5V, 500 μ s

TRIGGER IN socket allows the stimulator to be triggered from a positive going edge of an applied external signal, which should be TTL compatible.




5V

1.3 TECHNICAL SPECIFICATIONS

1.3.1 STIMULUS OUTPUT	
Current range	0 - 5 mA in 0.2 mA steps / 5 - 40 mA in 1 mA steps
Pulse width	100 - 1000 microseconds square wave (100 μ s step)
Maximum output voltage	50 V
Number of electrodes	9

1.3.2 STIMULUS PULSE RATES	
Pulse or group rate	0.1 to 8.0 per second, in 0.1 steps
Burst selections	Number of pulses - 1 to 9
Interstimulus time	1 - 10 ms in 0.1 ms steps



1.3.3 TRIGGER FACILITIES	
Manual	Push button mounted on front panel
Auto	
External	TTL compatible positive edge logic signal
Trigger output	Positive TTL compatible signal 500 μ s wide

1.3.4 OTHER	
Mechanical dimensions	250x263x85 mm
Weight	2 kg
Power	100 - 240 V @ 50-60 Hz
Rating	24.4-31.7 VA
Startup time	30 sec


1.4 MAINS CONNECTION

This device is shipped complete with mains lead if the mains plug has to be removed because of incompatibility with the local situation; when connecting another type of mains plug, the connections are:

Line - black
 Neutral - yellow
 Ground - green/yellow

A mains plug must be constructed so that if the lead is pulled out, the ground wire cannot come in contact with the line or neutral pins.

Fuse that is to be replaced must be of correct dimensions and be rated at 2 Amps.



Fuses that should be used when replacing are: 2 x High breaking capacity fuse 2 A (250)

Stimulator power supply module is of medical type, and can be used for mains voltages between 110 and 220 V (50 to 60 Hz) without additional adjustments (automatic voltage control)

1.5 ENVIRONMENTAL CONDITIONS

Transport and storage temperature: -19°C to 60°C

Transport and storage relative humidity: 10% to 80% (non-condensing)

Please Note: Following transportation and if storage temperature is outside the operating temperature range, allow the unit to acclimate for 3 hours prior to use.

Operating temperature: 5°C to 30°C

Operating relative humidity: 30% to 70% (non-condensing)

Atmospheric pressure: 50kPa to 106kPa

Maximum altitude: 3000 m (height above sea level)

The device does not deliver any heat to the environment.

1.6 OUTPUT IMPEDANCES AND CURRENT

VOLTAGE (V)	IMPEDANCE (OHM)	CURRENT (mA)
50	500	40
50	1000	40
50	2000	25
50	3000	16
50	4000	12

1.7 ELECTRODES INTENDED FOR USE

There are several disposable neurological electrodes intended for use with this stimulator for a specific purpose. Only electrodes terminated with Safety DIN42802 connectors should be used with the stimulator. Electrodes are not included in SGM-STIM-2 package. Electrodes must be verified for use in IONM. Any electrodes used, different from the ones stated above, are excluded as they may cause unforeseeable risk.




2. INSTRUCTIONS FOR USE

2.1 PRECAUTIONS WHEN USING THE SGM-STIM-2 STIMULATOR









The following remarks apply to the use of stimulators designed for percutaneous use: To avoid any risk of establishing an unintentional path to earth, connections from the stimulator to the patient have to be isolated from the ground. This means that a ground connection must not be made to either the terminal on the stimulator or to stimulating electrode. It is also undesirable to place a ground electrode on the patient close to one of the stimulating electrodes. A patient with an implanted electrode device should not be subjected to electrical stimulation unless specialist medical opinion has been obtained first.

Explanations of symbols used on SGM-STIM-2

Front Panel

	ATTENTION, CONSULT ACCOMPANYING DOCUMENTS
	CAUTION
	TYPE BF APPLIED PART

Back Panel

	TRIGGER INPUT
	TRIGGER OUTPUT
	POSITIVE GOING PULSE
	SQUARE WAVE PULSE
	OFF (POWER)
	ON (POWER)
	EQUIPOTENTIALITY
	DATE OF MANUFACTURE

2.2 POWERING THE STIMULATOR

Before powering the stimulator, toggle switch on the front panel has to be set in HIGH VOLTAGE OFF position. Connect the stimulator to the mains supply with an appropriate cable. Turn the stimulator on with the main power switch on the rear panel. Wait approximately 30 seconds to allow internal microprocessor and high-voltage generator to be operable. After switching POWER-ON and a short period of booting, the following factory preset menu appears on the display.

Select mode:

- multipulse strip
- multipulse handheld
- Penfield strip
- Penfield handheld

2.3 DEFINITION OF STIMULUS PROPERTIES

Stimulator output waveform is a rectangular pulse train with adjustable parameters as defined in the section "1.3 TECHNICAL SPECIFICATIONS". Using display and navigation keyboard, select output current amplitude, number of output pulses in a group, stimulus ratio and inter-stimulus interval.

2.4 OPERATING MODE SETUP

"Trigger:" display entry is used to select trigger mode: manual, external or auto. In "auto" mode, pulse train repetition rate can be adjusted between 0.1 and 8.0 Hz.

2.5 CONNECTING THE PATIENT

Connect stimulator output connector to the corresponding strip electrode and neutral electrode montage on the patient head. Use approved electrodes with appropriate cabling. Select active electrode using "output:" entry on the display. If the distance between patient and stimulator is exceeding the length of the electrodes, please use Output extender. See - Appendix 1. Stimulating electrodes must be selected in accordance with the specification in the item 2.8.

2.6 STIMULATION

Set the toggle switch HIGH VOLTAGE to ON position. In manual operating mode, push PULSE-START-STOP button to send single stimulus to the patient. In automatic operating mode, push PULSE-START-STOP button to start or stop generation of output stimuli.

2.7 END OF STIMULATION

1. Stop the generation of output pulses if in AUTO mode of operation, using PULSE-START-STOP button.
2. Set toggle switch HIGH VOLTAGE to OFF position.
3. Deselect active output electrodes by putting them in "off" mode.
4. Disconnect the electrodes from the patient.
5. Turn off the stimulator using the main switch.
6. Disconnect the mains cable connector from the power outlet.

2.8 STIMULATING ELECTRODES

Stimulating electrodes are all standard electrodes for intraoperative monitoring with impedance < 1 kOhm and connector DIN 42802.

2.9 OUTPUT LIMITS

The SGM-STIM-2 stimulator is a constant current type stimulating device, so it is relatively insensitive to the variation of the load impedance, as long as this impedance is within reasonable limits. These limits are defined with the maximum allowed output voltage that can be applied to the patient. Maximum allowed output voltage of the stimulator determines maximum current which can be delivered into load impedance, according to the following relation:

$$I_{out_max} = 210V / Z_{out}$$

OUTPUT IMPEDANCE	MAXIMAL OUTPUT CURRENT
0,5 kOhm	40 mA limited by setup
1 kOhm	40 mA limited by setup
2 kOhm	25 mA limited by maximum output voltage
3 kOhm	16 mA limited by maximum output voltage
4 kOhm	12 mA limited by maximum output voltage

Maximum output voltage that stimulator can deliver can be factory preset according to customer requirements and FDA or IEC regulations.

2.10 DEVICE INSPECTION AND MAINTENANCE

Routine check of the stimulator is performed after any powering of the device. Factory preset values as defined in Section 2.2 must appear on the display. Accuracy of operating data can be checked by qualified personnel using 1 kOhm load resistance and appropriate equipment for measuring intensity and duration of output current pulses. The values of pulse durations, pulse repetition frequencies and amplitudes, as indicated on the equipment, shall not deviate by more than 30%, when measured with an error not exceeding +/-10%.

Device has no wearable parts that are to be replaced regularly. In case of appearance of the DC component on the output of stimulator, immediately switch off the preset pair of electrodes on the output selector. Preventive service is obligatory at least once every 3 years. SGM d.o.o. is the only one authorized for service and calibration of device. No modifications of this equipment are allowed.


2.11 CLEANING

Before cleaning, disconnect the device from mains outlet. Cleaning of the system should be kept to the minimum. If cleaning is needed, it should be done very carefully due to the sensitivity of the system. Use non-alcohol based disinfectant wipes. Do not use: alcohol, thinners or any other solvent to clean the system. Do not use abrasive, corrosive cleaning agents or solvents. Do not sterilize any components.

3. ACCESSORIES

Output Extender:

Output extender is designed for easier operation with the patient, to keep the medical stimulating equipment in appropriate distance from the patient. It is designed to easily identify and connect the proper output contact from the stimulator and to transfer it on the mobile, hand-held panel. The cable terminates with DIN 42802 touch proof connector. Each connector is marked according to the labeling on the front panel of the stimulator and output extender panel. On the back panel of the output extender is the clip for an easier stabilization of the panel in the place nearest to the patients.



Pedal:

Foot pedal has the same function as the front panel mounted yellow button, marked as PULSE – START – STOP. 1.5 mm termination of the foot pedal cable should be connected to the back panel of the device, in the socket marked PEDAL.

4. EMC TABLE

THE DEVICE IS INTENDED FOR USE IN THE ELECTROMAGNETIC ENVIRONMENT SPECIFIED BELOW. THE CUSTOMER OR THE USER OF THE DEVICE SHOULD ASSURE THAT IT IS USED IN SUCH AN ENVIRONMENT.			
Emissions test	Compliance	Electromagnetic environment - guidance	
Mains Terminal Spurious Voltage EN 55011-1 Clause 8.2	Class B	The lower limit is valid for the frequency boundaries	
Radiated Field EN 55011-1 Clause 8.2	Class A	The device is suitable for use in all establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.	
Harmonic Current Emissions 61000-3-2	Class A		
Voltage Variation Reduction and Flickering Reduction in the Mains IEC 61000-3-3	Complies		
IMMUNITY TEST	IEC 60601 TEST LEVEL	COMPLIANCE LEVEL	ELECTROMAGNETIC ENVIRONMENT - GUIDANCE
Electrostatic discharge (ESD) IEC 61000-4-2	±4 kV contact ±4 kV air	PERFORMANCE CRITERION B ±4 kV contact ±4 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%. Contact discharge ±4 kV could cause display reset.
Field EN 61000-4-3	80 MHz – 2.5 GHz	PERFORMANCE CRITERION A	Mains power quality should be that of a typical commercial or hospital environment.
Electrical fast transient/ burst IEC 61000-4-4	±0,5 kV for power supply lines ±0,25 kV for input/ output Extender cable	PERFORMANCE CRITERION A ±0,5 kV for power supply lines ±0,25 kV for input/ output Extender cable	Mains power quality should be that of a typical commercial or hospital environment. Temporary loss of charging function may occur, which is self-recoverable.
Surge IEC 61000-4-5	±0,5 kV L conductor – PE conductor ±0,5 kV N conductor – PE conductor	PERFORMANCE CRITERION A ±0,5 kV L conductor – PE conductor ±0,5 kV N conductor – PE conductor	Mains power quality should be that of a typical commercial or hospital environment.

Conducted Disturbances Induced by RF Fields IEC 61000-4-6	Conductors L, N and PE 150 kHz-80 MHz Output Extender cable 150 kHz-80 MHz	PERFORMANCE CRITERION A Conductors L, N and PE 150 kHz-80 MHz Output Extender cable 150 kHz-80 MHz	Mains power quality should be that of a typical commercial or hospital environment.
Power Frequency Magnetic Field IEC 61000-4-8	Frequency 50 Hz Field intensity 3	PERFORMANCE CRITERION A	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<0 % U _T for 0,5 cycle	<0% U _T for 0,5 cycle PERFORMANCE CRITERION A	Mains power quality should be that of a typical commercial or hospital environment.
	40% U _T for 5 cycle	40% U _T for 5 cycle PERFORMANCE CRITERION A	
	70% U _T for 25 cycle	70% U _T for 25 cycle PERFORMANCE CRITERION A	
	0% U _T for 250 cycle	0% U _T for 250 cycle PERFORMANCE CRITERION A	
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	1 Vrms 150 kHz to 80 MHz	[V1] 1 V	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = \left[\frac{3,5}{V_1} \right] \sqrt{P}$
Radiated RF IEC 61000-4-3	3V/m 80 MHz to 2,5 GHz	[E1] 3 V/m	$d = \left[\frac{3,5}{E_1} \right] \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = \left[\frac{7}{E_1} \right] \sqrt{P} \quad 800 \text{ MHz to } 2,5 \text{ GHz}$ where P is maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).

Neurostimulator SGM-STIM-2 in operation does not use any type of EMF-based device for remote control of its functionalities.

Metal grounded housing of the stimulator behaves as a Faraday cage, blocking any transmission of EM fields, so functionality of SGM-STIM-2 stimulator cannot be affected by external EMF.

TROUBLESHOOTING

To protect patients from any uncontrolled stimulating activity in case of a stimulator control computer dysfunctionality, SGM-STIM-2 has a protecting system which blocks the device high voltage output if control computer does not run its software.

ISSUE	POSSIBLE CAUSE	POSSIBLE SOLUTION
Neurostimulator does not start	Problem with control computer or power supply module	Return device to the producer for testing and servicing
Device does not deliver stimulus to the patient	Stimulator output is not correctly selected, HV switch is in OFF position	Use correct output setup procedure for SGM-STIM-2
Device does not deliver stimulus to the patient	Electrode is disconnected from the patient	Check electric connection path from stimulator output to the patient. Check implantation of stimulating electrode to the patient
Stimulator triggering with external IOM device does not work	External trigger impulse amplitude or duration are not within SGM-STIM-2 specification	Correct external trigger impulse parameters, if possible

Manufactured by:



Address: Grge Novaka 22a, 21000 Split, Croatia

Tel: ++385 (0) 21 779 916

Email: sales@sgm.hr

Web: www.sgm.hr